





## Products from Bandage Cellulose, Highly Bleached

### Description and application:

- Made of highly bleached cellulose fibres in the form of several creped single layers which are laid one above the other and felted with each other, TCF (totally chlorine-free bleached), meeting the requirements of DAB (monograph of highly bleached bandage cellulose - cellulosum ligni depuratum);
- Soft and abrasion-proof, can be used in various areas of medicine, in particular for wound care in combination with other wound care systems as an additional, cost-effective absorbent, protective and padding material; in the case of minor injuries as absorbant padding and for "minor disinfection".

---

## Bandage Cellulose, Pre-Packed

### Packaging:

- Delivered in printed flap bags, ready for sale in a pharmac

[read more...](#)

---

## Cellulose Pads

### Packaging:

- In special collapsible boxes;
- In printed flap bags

[read more...](#)

---

## **Bandage Cellulose Cuts**

### **Packaging:**

- 1,000 g original packs in a flap bag;
- 5 kg and 10 kg original packs in a cardboard box, in part as rolls wrapped in foil, ready for sale in a pharmacy

[read more...](#)