





Absorbent Compresses

- Made of different material layers;
- Various dimensions;
- Fitted with an absorbent core.
- The material layer which comes into contact with the wound surface displays good wound compatibility, i.e. no sticking to the wound.
- The compress rear side is coloured for simple identification and is furthermore water-resistant to prevent penetration of the wound secretion.
- Thanks to its good absorbency, the compress is the ideal choice for strongly secretory wounds.

Packaging:

- Non-sterile 10 pcs. in a pressure seal bag
20 or 25 pcs. in folding box
- sterile individually sealed

Size	Original Packaging	Article No.	PZN	Shipping Unit
10 x 10 cm				
Non-sterile	10 pcs.	132401	07354066	50 pack/ct.
	25 pcs.	132402	07354043	32 pack/ct.
Sterile	2 x 1 pcs.	132417	08828796	As required
	5 x 1 pcs.	132416	08828804	30 pack/ct.

Size	Original Packaging	Article No.	PZN	Shipping Unit
	25 x 1 pcs.	132415	08828810	6 pack/ct.
10 x 20 cm				
Non-sterile	10 pcs.	132405	07353782	50 pack/ct.
	25 pcs.	132406	07350186	12 pack/ct.
Sterile	2 x 1 pcs.	132418	08828827	As required
	5 x 1 pcs.	132419	08828879	30 pack/ct.
	25 x 1 pcs.	132420	08828939	6 pack/ct.
15 x 25 cm				
Non-steril	10 pcs.	132403	14293282	25 pack/ct.
	20 pcs.	132404	14293299	12 pack/ct.
Sterile	2 x 1 pcs.	132421	14293230	As required
	5 x 1 pcs.	132422	14293247	30 pack/ct.
	25 x 1 pcs.	132423	14293253	6 pack/ct.