



-



Dental Wadding Rolls

- Made of 100% cotton wool, Ph.Eur.;
- Important aid in dental surgery;
- Especially suitable as absorbent material to keep the treated area dry.

Packaging:

- In a practical, reclosable PE bag

Size	Original Packaging	Article No.	PZN	Shipping Unit
Size 1	Ø 8 mm, ca. 300 g	121210	07314492	10 pack/ct.
Size 2	Ø 10 mm, ca. 300 g	121211	07314500	10 pack/ct.
Size 3	Ø 12 mm, ca. 300 g	121212	07314517	10 pack/ct.