





## FES-TUPF®

### Properties and special features:

- Cut edges turned safely to the inside thanks to the special FESMED folding;
- Fast absorption of all kinds of wound fluid even around corners;
- Maximum absorbency if two pieces are used;
- Optimised, size-dependent absorbency if one piece is used;
- Tweezers can "reach into" the swab;
- Suitable for large-area skin disinfection;
- Highly flexible in use thanks to the large variety of sizes

Size	Original Packaging	Article No.	PZN	Shipping Unit
<b>No. 3, lentil-size</b>				
Non-sterile	500 pcs.	162150	04050621	As required
Sterile	5 x 10 pcs.	162139	14266345	20 pack/ct.
	10 x 10 pcs.	162140	14266411	10 pack/ct.
<b>No. 4, pea-size</b>				
Non-sterile	500 pcs.	162151	04050638	As required
Sterile	5 x 10 pcs.	162141	14266428	20 pack/box
	10 x 10 pcs.	162142	14266434	10 pack/box
<b>No. 6, bean-size</b>				
Non-sterile	100 pcs.	162159	06310747	As required
	400 pcs.	162152	04050644	As required
Sterile	5 x 5 pcs.	162143	14266440	15 pack/box
	10 x 5 pcs.	162144	14266457	10 pack/box
<b>No. 9, hazelnut-size</b>				
Non-sterile	100 pcs.	162160	06310753	As required
	300 pcs.	162153	04050561	As required
Sterile	5 x 5 pcs.	162180	10415179	15 pack/box
	100 x 5 pcs.	162181	10415185	12 pack/ct.

Size	Original Packaging	Article No.	PZN	Shipping Unit
	5 x 25 pcs.	162145	14266463	5 pack/box
	2 x 75 pcs.	162146	14266486	5 pack/box
	10 x 75 pcs.	162147	14266500	16 pack/ct.
No. 12, walnut-size				
Non-sterile	100 pcs.	162161	06310776	20 bag/ct.
	500 pcs.	162154	04050650	10 bag/ct.
Sterile	5 x 1 pcs.	162133	11158589	10 pack/box
	20 x 1 pcs.	162175	10415216	5 pack/ct.
	100 x 1 pcs.	162176	10415222	12 pack/ct.
	5 x 2 pcs.	162174	10305031	10 pack/box
	50 x 2 pcs.	162172	00261396	16 pack/ct.
	25 x 3 a 5 pcs.	162173	00261404	8 pack/ct.
	5 x 6 pcs.	162134	11158595	5 pack/box
	20 x 6 pcs.	162177	10415239	12 pack/ct.
	5 x 15 pcs.	162135	11158603	3 pack/box
	15 x 15 pcs.	162178	10415251	12 pack/ct.
	2 x 50 pcs.	162169	14266598	4 pack/box
	10 x 50 pcs.	162170	14266606	12 pack/ct
No. 15, plum-size				
Non-sterile	100 pcs.	162162	06310782	20 bag/ct.
	500 pcs.	162156	04094512	10 bag/ct.
Sterile	5 x 1 pcs.	162184	11158661	10 pack/box
	20 x 1 pcs.	162185	10415268	5 pack/box
	100 x 1 pcs.	162186	10415274	12 pack/ct.
	5 x 2 pcs.	162182	11158649	10 pack/box
	50 x 2 pcs.	162183	11158655	16 pack/ct.
	5 x 5 pcs.	162187	11158626	4 pack/box
	20 x 5 pcs.	162188	10415280	12 pack/ct.
	5 x 6 pcs.	162189	11158632	4 pack/box
	20 x 6 pcs.	162190	10415297	12 pack/ct.
	5 x 12 pcs.	162192	10415305	As required
	25 x 12 pcs.	162193	10415311	8 pack/ct.
	2 x 30 pcs.	162167	14266612	4 pack/box
	10 x 30 pcs.	162168	14266629	12 pack/ct.