

•



•



•



## Overview of Gauze Swabs

**Our range of supply includes the following variants:**

- FES-TUPF®- triangular hat-style form based on the FESMED technology;
- Ball-shaped gauze swabs, cut edges held inside tightly by way of a rubber ring;
- Ball-shaped gauze swabs, without a rubber ring

**Description and application:**

- 100% cotton, aseptic gauze, DIN EN 14079, Ph. Eur., type 20;
- to halt bleeding and for surgical cleaning, as well as for large-area skin disinfection.

**Packaging:**

- Non-sterile: Packed in pressure seal bags;
- Sterile: Sealed in special-paper bags with sterilisation indicator, packed in a collapsible box

**Sizes:**

The various size designations found on the market are compared below to facilitate orientation. The number of the FESMED designation is used as the reference size.

Number	Description	FES-TUPF® Triangular, hat-style form	Ball-shaped gauze swabs, without rubber ring Ø	Ball-shaped gauze swabs, with rubber ring Ø
3	Lentil-size	Approx. 1.0 cm		
4	Pea-size	Approx. 1.3 cm		
6	Bean-size	Approx. 1,5 cm		
9	Hazelnut-size	Approx. 2,0 cm		Approx. 2,0 cm

Number	Description	FES-TUPF® Triangular, hat-style form Side length	Ball-shaped gauze swabs, without rubber ring Ø	Ball-shaped gauze swabs, with rubber ring Ø
12	Walnut-size	Approx. 3,5 cm		Approx. 2,5 cm
15x15	Walnut-size		Approx. 2,5 cm	
20x20	Plum-size		Approx. 3,0 cm	
15	Plum-size	Approx. 5,0 cm		Approx. 3,5 cm
15b	Egg-size		Approx. 4,0 cm	
30	Extra-large		Approx. 6,0 cm	
40	Extra-large		Approx. 7,0 cm	