



-



Tamponades

- Selvedge cellulose gauze tape, glossing, rolled up to a width of 2 cm, folded from a width of 6 cm;
- To be used for padding, tamponading and inserting strips into small body cavities;
- In the oral cavity up to the throat, in the auditory canal up to the eardrum or into the nasal cavity;
- For transient use.

Packaging:

- Non-sterile: in a special collapsible box;
- Sterile - sealed in special-paper bags with sterilisation indicator, packed in collapsible box

Size	Original Packaging	Article No	PZN	Shipping Unit
1 cm x 2,5 m wound double				
Non-sterile	10 pcs.	171210	07402569	As required
Sterile	5 x 1 pcs.	171222	00263633	10 pack/box
	50 x 1 pcs.	171220	00263515	16 pack/ct.
1 cm x 5 m				
Non-sterile	10 pcs.	171211	07402575	As required
Sterile	5 x 1 pcs.	171225	00264940	10 pack/box
	50 x 1 pcs.	171223	00263656	16 pack/ct.
2 cm x 5 m				
Non-sterile	10 pcs.	171212	07402581	As required
Sterile	5 x 1 pcs.	171227	00263952	10 pack/box
	25 x 1 pcs.	171226	00263946	16 pack/ct.
6 cm x 5 m				
Sterile	2 x 1 pc.	171218	11158678	10 pack/box
8 cm x 5 m				
Sterile	2 x 1 pcs.	171219	14040977	10 pack/box